

TYPE OF ANAESTHESIA

Submucosa papillary area	Para-apical by vestibular path Upper and lower jaw	Para-apical by palatine path Upper jaw	Para-apical by lingual path Lower jaw	Intra-osseous Intra-septal	Intra-ligamentous	Intra-pulpal	Regional Anaesthesia Inferior Alveolar Nerve Block
--------------------------	--	--	---------------------------------------	----------------------------	-------------------	--------------	--

TYPE OF NEEDLE

●	30 G X-SHORT 0,30 x 12 mm	●●●	●	●●	●	●●●	
	30 G SHORT 0,30 x 16 mm	●	●●●	●●●	●●●	●	●
	30 G SHORT 0,30 x 10 mm		●●●	●●●	●●●		●●●
	30 G SHORT 0,30 x 23 mm			●●	●●	●	●●
	30 G SHORT 0,30 x 25 mm		●●	●●		●	
●	27 G X-SHORT 0,40 x 8 mm	●			●●●	●●●	
	27 G SHORT 0,40 x 16 mm		●●	●●	●	●	
	27 G SHORT 0,40 x 21 mm		●●	●●	●		●
	27 G SHORT 0,40 x 23 mm		●	●	●		●
	27 G SHORT 0,40 x 25 mm		●	●	●		
	27 G SHORT 0,40 x 32 mm		●	●			
	27 G-LONG 0,40 x 35 mm						●●
●	25 G-LONG 0,50 x 35 mm						●●

●●● Recommended needle ●● Useable needle ● Possible use