



by **SOFIC**

injection needle

AVAILABLE SIZES

TYPE OF ANAESTHESIA

TYPE OF NEEDLE		Submucosa Papillary area	Peri-apical vestibular path Maxil. & Mandib.	Maxilar peri-apical palatine path Maxilla	Peri-apical lingual path Mandible	Intra-osseous Intra-septal	Intra-ligamentous	Intra-pulpal	Regional Anaesthesia Inferior Alveolar Nerve Block
		▼	▼	▼	▼	▼	▼	▼	▼
30 G X-SHORT	0,30 x 12 mm	●●●	●	●●		●	●●●		
30 G SHORT	0,30 x 16 mm	●	●●●	●●●	●●●		●	●	
30 G SHORT	0,30 x 21 mm		●●●	●●●	●●●			●●●	
30 G SHORT	0,30 x 23 mm		●●	●●	●			●●	
30 G SHORT	0,30 x 25 mm		●●	●●				●	
27 G X-SHORT	0,40 x 8 mm	●				●●●	●●●		
27 G SHORT	0,40 x 16 mm		●●	●●	●	●	●		
27 G SHORT	0,40 x 21 mm		●●	●●	●			●	
27 G SHORT	0,40 x 23 mm		●	●	●			●	
27 G SHORT	0,40 x 25 mm		●	●	●				
27 G-LONG	0,40 x 32 mm		●	●					●●
27 G-LONG	0,40 x 35 mm								●●

●●● Recommended needle ●● Useable needle ● Possible

Private label available

EC CERTIFICATE
N° 01085



ISO 14001
EMS 39265



ISO 9001
FM 32741

ISO 13485
MD 80441
FM 86084

SOFIC ● 3, rue Jean-Jacques Rousseau ● BP 282 Aussillon ● 81209 Mazamet Cedex - France
Tel : + 33 (0) 5 63 97 71 65 ● Fax : +33 (0) 5 63 97 71 66 ● Email = info@sofic.com

SA au capital de 912 003 Euros - RC Castres B 716.720.248 - Siret 716.720.248.00039 - Code Naf 331 B - N° TVA Intracommunautaire FR 71 716.720.248

www.sofic.com