

## ACCESSORIES DESIGNED FOR ALL YOUR INJECTION REQUIREMENTS

### PLASTIC HANDLE



Black ring  
Ref: PPAU

The SOFICONCEPT XL single use, safety syringe comes with a black ring handle as standard, one to a box of 100 units. Made from polycarbonate, this is non-reflective and warm to the touch and the design allows for both passive and active aspiration. A great all rounder that lasts for a minimum of one hundred autoclave cycles and can also be bought separately from the syringe system. However, there are now three other choices available to the dentist: the METAL RING, T-BAR and SOFINJECT.

### METAL HANDLES



Metal ring\*  
Ref: PPAT2

The METAL RING and T-BAR handles are designed for those who prefer something a little heavier than the black ring handle. They are made from stainless steel and polycarbonate.

METAL RING handle is suitable for passive and active aspiration.

T-BAR handle is suitable for passive aspiration only.



T BAR\*  
Ref: PPAT1

The SOFINJECT handle allows the practitioner to carry out intraligamentary and intra-osseous injections whilst still enjoying the safety that SOFICONCEPT XL provides. Its ergonomic design and smooth action allow the practitioner to deliver a slow injection into areas where strong pressure is often needed. The forward movement of the piston releases 0.2 ml of anaesthetic each time the handle is moved so the practitioner knows exactly how much anaesthetic has been given.



SOFINJECT\*  
Ref: PPADN

All the handles are designed to work exclusively with the SOFICONCEPT XL syringe and when attached correctly it takes 40 kilos of pressure to part the handle from the barrel. The handles are all fully autoclavable.

\*Optional accessory