



TECHNICAL SPECIFICATIONS

● ISO 9626 stainless tube



◀◀ ● Triple bevel

● Silicone coating

● Threaded hub : metric or inch thread

● Transparent sheath/cap



◀◀ ● Bevel indicator\*

● Ethylene oxide sterilization

\*only available on customer's demand



## AVAILABLE SIZES

### TYPE OF ANAESTHESIA

TYPE OF NEEDLE	Needle Size	Submucosa Papillary area	Peri-apical vestibulary path Maxil. & Mandib.	Maxilar peri-apical palatine path Maxillia	Peri-apical lingual path Mandible	Intra-osseous Intra-septal	Intra-ligamentous	Intra-pulpal	Regional Anaesthesia Inferior Alveolar Nerve Block
		●●●	●●	●●	●●	●●	●●	●●	●●
●	<b>30 G X-SHORT</b> 0,30 x 8 mm	●●●	●			●	●●●		
	<b>30 G X-SHORT</b> 0,30 x 12 mm	●●●	●	●	●		●●		
	<b>30 G SHORT</b> 0,30 x 16 mm	●	●●●	●●●	●●●			●	
	<b>30 G SHORT</b> 0,30 x 21 mm		●●●	●●●	●●●			●●●	
	<b>30 G SHORT</b> 0,30 x 23 mm		●●	●●	●			●●	
	<b>30 G SHORT</b> 0,30 x 25 mm		●●	●●	●			●	
●	<b>27 G X-SHORT</b> 0,40 x 8 mm	●				●●●	●●●		
	<b>27 G SHORT</b> 0,40 x 16 mm		●	●	●	●			
	<b>27 G SHORT</b> 0,40 x 21 mm		●●	●●	●●			●	
	<b>27 G SHORT</b> 0,40 x 23 mm		●●	●●	●			●	
	<b>27 G SHORT</b> 0,40 x 25 mm		●●	●●	●				
	<b>27 G-LONG</b> 0,40 x 32 mm		●	●					●
	<b>27 G-LONG</b> 0,40 x 35 mm								●●●
	<b>27 G-LONG</b> 0,40 x 42 mm								●
●	<b>25 G-LONG</b> 0,50 x 35 mm								●●●
	<b>25 G-LONG</b> 0,50 x 42 mm								●●

●●● Recommended needle    ●● Useable needle    ● Possible

● Private label available

**EC CERTIFICATE**  
N° 01085



**ISO 14001**  
EMS 39265



**ISO 9001**  
FM 32741  
**ISO 13485**  
MD 80441  
FM 86084

SOFIC ● 3, rue Jean-Jacques Rousseau ● BP 282 Aussillon ● 81209 Mazamet Cedex - France  
Tel : + 33 (0) 5 63 97 71 65 ● Fax : +33 (0) 5 63 97 71 66 ● Email = info@sofic.com

SA au capital de 912 003 Euros - RC Castres B 716.720.248 - Siret 716.720.248.00039 - Code Naf 331 B - N° TVA Intracommunautaire FR 71 716.720.248

[www.sofic.com](http://www.sofic.com)